

recountmethylation cheatsheet

TERMINOLOGY

Download DNAm compilations

DNA methylation (DNAm) terms

CpG locus	: Cytosine-guanine dinucleotide pair
DNA methylation (DNAm)	: DNA with covalently bound methyl (CH_3) groups; here specify cytosine-bound methyl groups in a CpG locus.
CG Island	: DNA regions enriched for cytosine, guanine, and CpG loci.

BeadArray terms

CpG probe	: Class of DNAm array technology using bead-bounded probes to quantify DNAm at specific CpG loci.
BeadArray	: Type of microarray using BeadArray probes, manufactured by Illumina, to quantify DNAm.
HM27k	: Type of DNAm array, introduced around 2005, targeted roughly 27,000 CpG loci in humans.
HM450k	: Type of DNAm array, introduced around 2011, targeting roughly 450,000 CpG loci in humans.
EPIC	: Type of DNAm array, introduced around 2015, targeting roughly 855,000 CpG loci in humans. Shares 93% (about 453,000) of probes with the HM450k platform.

Data object classes

The following command line options show how to download database files from the server (<http://www.recount.bio/data> (<http://www.recount.bio/data>)), or you may try to right-click -> download from within your browser. Note, it may help to increase your timeout period for long downloads.

- Use `recountmethylation`: e.g. `getdb_h5se_rg("hm450k")` to download the HM450k `RGChannelSet` data as an `HDF5-SummarizedExperiment` object (see `?getdb` for more info).
- Use `wget`: From command line, enter `wget -r <filepath>`, replacing `<filepath>` with an address from <https://www.recount.bio/data>. Note the `-r` is needed for `h5se` objects, which are directories.

DelayedArray operations and pipes

The following operations make use of `DelayedArray` caching.

- Perform summary operations with `dim()` / `nrow()` / `ncol()`
- Rapidly update the sample metadata in an `h5se` object with `quickResaveHDF5SummarizedExperiment()`.
- Automatically pipe data chunks between `DelayedArray` objects, e.g.

```
library(HDF5Array)
h5 <- loadHDF5SummarizedExperiment(h5se.path)
bval <- getBeta(h5) # this is delayed
saveHDF5Array(bval, file = "db.h5") # executes with chunking
```

HDF5 / h5	: Hierarchical database format 5, a type of database syntax implementing compression and chunking.
-----------	--

RECAST DATASET OBJECTS

se	: Short for <code>SummarizedExperiment</code> , a multifaceted object class in R/Bioconductor containing slots for assay measurements and metadata for the platform, samples, and experiment.
----	---

If operations on datasets throw errors due to their class, you may attempt to recast and rerun them without the `DelayedArray` backend.

Recast an `RGChannelSet`

```
library(minfi)
rg <- loadHDF5SummarizedExperiment(rg.path)
green.matrix <- as.matrix(getGreen(rg))
red.matrix <- as.matrix(getRed(rg))
anno <- annotation(rg)
metadata <- DataFrame(pData(rg))
rg.new <- RGChannelSet(Green = green.matrix,
                      Red = red.matrix,
                      annotation = anno,
                      colData = metadata)
```

Recast a `GenomicRatioSet`

```
library(minfi)
gr <- loadHDF5SummarizedExperiment(gr.path)
bval.matrix <- as.matrix(getBeta(gr))
anno <- annotation(gr)
metadata <- DataFrame(pData(gr))
ranges <- granges(gr)
gr.new <- GenomicRatioSet(gr = ranges,
                        Beta = bval.matrix,
                        annotation = anno,
                        colData = metadata)
```

h5se	: Short for <code>HDF5-SummarizedExperiment</code> , a hybrid object class that uses <code>DelayedArray</code> for caching.
RGChannelSet / rg	: Type of <code>se</code> object containing dual color channel data for red and green channels on Illumina's BeadArray DNAm array platforms.
GenomicMethylSet / gm	: Type of <code>se</code> object containing the methylated (a.k.a. 'M') and unmethylated (a.k.a. 'U') signals calculated from dual channel intensity data.
GenomicRatioSet / gr	: Type of <code>se</code> object containing the DNAm fractions (a.k.a. Beta-values) and/or logit2-transformed fractions (a.k.a. M-values) calculated from the M and U fractions.

DECISION TREE

Once you've selected a DNAm array platform (e.g. HM450k or EPIC), you can use the following decision tree to determine which database compilation object type to download.

